

Sexing bovine preimplantation embryos by PCR

Yao Xiao, Froylan Sosa and Peter J Hansen

Department of Animal Sciences, University of Florida

Introduction

This protocol for using PCR to sex bovine embryos are based on Y-specific and autosomal targets described by Park et al, Theriogenology 2001; 55: 1843-1853 but with updates according to current bovine genomic database in NCBI as of 2019. Note that assembly of the bovine Y chromosome is still in progress.

Primers

Y specific primers

The Y-specific region targeted for PCR is BOV97M located at CM001061.2 (GenBank ID):

```
gatcactatacacaccactctcatcctacctaataagattccagctgataattgcaactatactccaaagtttgaggaaatgtgaaagaaaatc  
attcagctcatctgtttcagagacagaatttgtgtagcatagcctttttatcacaagatc
```

Forward: 5'-gatcactatacacaccact-3'

Reverse: 5'-aaggctatgctacacaaattctg-3'

size of amplicon = 143bp

Targeted autosome: Chromosome 19

Autosomal primers (5' → 3'), size of amplicon 217bp*:

Forward: 5'tggaagcaaagaaccccgct-3'

Reverse: 5'tcgtgagaaaccgcaccctg-3'

Size of amplicon = 217bp (other bands at 1623 and 3029bp).

Embryo Collection

1. Wash blastocysts with DEPC-treated Dulbecco's phosphate buffered saline containing 0.2% (w/v) polyvinylpyrrolidone (DPBS-PVP) three times.
2. Remove zona with Acid Tyrode's (Sigma T1788).
3. Repeat step 1.
4. Move a blastocyst with minimum amount solution (<0.5 µl) into a PCR tube with 1 µl of resuspension buffer from the CellsDirect kit 11753 (ThermoFisher Scientific).

- Snap freeze and store at -80°C.

Note: If you work with hatched blastocysts (blastocysts older than 7 days-old), we strongly recommend exposing embryos to drops of 50 µL Dissociation Reagent (TrypLE™ Select 10x, Gibco, REF: A12177-01) to ensure full cell digestion. Exposure to Acid Tyrode's is not required when using the Dissociation Reagent.

Embryo lysis and RNase treatment

- Add 0.5 µl of lysis enhancer from the CellsDirect kit 11753 (ThermoFisher Scientific) into each tube containing an embryo under a stereoscope.
- Digest at 70°C for 20 min in a thermocycler. At 10 min, tap the tubes to mix and then continue digestion.
- Make sure embryos are fully dissolved under a stereoscope. If not, digest for 10 min more at 70°C.
- Add 1 µl of 100 µg/mL RNase A (19101, Qiagen; diluted with nuclease free water). Incubate at 37°C for 0.5 h.

PCR

- Prepare 10mM dNTP:
Add 10 µl of dATP, dGTP, dCTP, dTTP (10297-018, Invitrogen) into 60 µl nuclease free water.

- Prepare mastermix as follow:


Items	Volume (µl)
5X Green GoTap Flexi Buffer (M891A, Promega)	4
25 mM MgCl ₂ (A351H, Promega)	4
10 mM dNTP	1
10 µM sex primers (forward + reverse)	0.4
GoTaq hot start polymerase (M500B, Promega)	0.2
Nuclease-free water	6.6
TOTAL	16.2

- Add 16.2 µl of mastermix into each tube containing a single embryo. The final volume will be ~ 19 µl.

4. As a positive control, add 16.2 μl of mastermix to a tube containing 2.8 μl of 10 ng/ μl genomic DNA isolated from blood of male and female.
5. As a negative control, add 16.2 μl of mastermix to 2.8 μl of nuclease free water.
6. 1st round PCR for sex primers: 95°C, 5 min; (95°C, 15s; 58°C, 15s; 72°C, 15s) X 20; 72°C, 10 min.
7. Add 1 μl of 4 μM of autosomal primer into each tube, and perform a 2nd round of PCR: 95°C, 5 min; (95°C, 15s; 58°C, 15s; 72°C, 15s) X 17; 72°C, 10 min.

Gel electrophoresis

1. Prepare a 2% (w/v) agarose gel: add 2 g of low EEO agarose (BP160, ThermoFisher) into 100 mL of TAE buffer 1X (diluted from 10X, T9650, Sigma). Heat in a microwave until it is clear.
2. Pour the gel with a comb, allow to solidify, and then remove the comb.
3. Submerge the gel beneath the TAE buffer.
4. Place 2 μl drops of Diamond Nucleic Acid Dye (H1181, Promega) onto a paraffin strip and mix with 8 μl of PCR products or PCR controls. For ladder: Place 2 μl drops of Diamond Nucleic Acid Dye onto a paraffin strip and mix with 10 μl DNA ladder.
5. Load the samples into the wells of the gel.
6. Run at 100 V, 160 mA for 40-60 min.
7. Take images using a photoimager.



Created August 6 2019 Modified October 1 2019